● ESPAÑOL

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
Use, Care, & Maintenance

LIMITED WARRANTY

American Diagnostic Corporation (ADC®) warrants its products against defects in materials and workmanship under normal use and service as follows:

1. Warranty service extends to the original retail purchaser only and commences with the date of delivery.
2. Your DIAGNOSTIX™ aneroid is warranted for one year from date of purchase (all parts).
3. The manometer is warranted to remain accurate to +/-3mmHg (or the prevailing standard) over its full range when compared to a reference standard for life.

What is Covered: Replacement of parts, and labor.

What is Not Covered: Transportation charges to and from ADC®. Damages caused by abuse, misuse, accident, or negligence. Incidental, special, or consequential damages. Some states do not allow the exclusion or limitation of incidental, special, or consequential damages, so this limitation may not apply to you.

To Obtain Warranty Service: Send item(s) postage paid to ADC®, Attn: Repair Dept., 55 Commerce Dr., Hauppauge, NY 11788. Please include your name and address, phone no., proof of purchase, and a brief note explaining the problem.

Implied Warranty: Any implied warranty shall be limited in duration to the terms of this warranty and in no case beyond the original selling price (except where prohibited by law).

This warranty gives you specific legal rights and you may have other rights which vary from state to state.

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To Register Your Product, visit us at

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www.adctoday.com

y siga los acoplamientos

FOR QUESTIONS, COMMENTS, OR SUGGESTIONS
CALL TOLL FREE:
1-800-ADC-2670

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www.adctoday.com
info@adctoday.com

ENGLISH

Thank you for choosing an ADC® DIAGNOSTIX™ Aneroid Sphygmomanometer. Your new ADC® DIAGNOSTIX™ aneroid is designed to provide unrivaled durability, and unparalleled comfort. Every component has been carefully designed to maximize performance.

Device Description and Intended Use

ADC® aneroid sphygmomanometers are used by professional healthcare providers and individuals trained in auscultatory blood pressure technique to determine systolic and diastolic blood pressure in humans and animals.


General Warnings

A warning statement in this manual identifies a condition or practice which, if not corrected or discontinued immediately could lead to patient injury, illness, or death.

- **WARNING:** Do not allow a blood pressure cuff to remain on patient for more than 10 minutes when inflated above 10 mmHg. This may cause patient distress, disturb blood circulation, and contribute to the injury of peripheral nerves.
- **WARNING:** If luer lock connectors are used in the construction of tubing, there is a possibility that they might be inadvertently connected to intra-vascular fluid systems, allowing air to be pumped into a blood vessel. Immediately consult a physician if this occurs.
- **WARNING:** Safety and effectiveness with neonate cuff sizes 1 through 5 is not established.
- **CAUTION:** For all blood pressure systems that can be wall mounted, ensure that the unit is securely mounted prior to use to avoid damage to the instrument and potential patient injury.
- **CAUTION:** Federal law restricts this device to sale by or on the order of a physician or licensed healthcare practitioner.
- **CAUTION:** To obtain the greatest accuracy from your blood pressure instrument, it is recommended that the instrument be used within a temperature range of 50°F (10°C) to 104°F (40°C), with a relative humidity range of 15%-90% (non-condensing).
- **CAUTION:** Unless otherwise indicated, this product contains natural rubber latex, which may cause allergic reactions.

Operation of Diagnostix™ 703, System 3, Pediatric, General Practice, System 4, and System 5™ Multicuff BP instruments

Your DIAGNOSTIX™ 703 consists of a rugged, palm style aneroid manometer with integral bulb and trigger air release valve (#804N), ADCUFF™ calibrated cuff with Size Guide™ marking system, bladder, luer slip connector (attached to bladder tubing), zippered carrying case, and instruction booklet. ADC®'s Multicuff Systems (System 3, 4, 5, GPK, and PMK) contain the #804 manometer and a variety of cuffs with single tube bladders, each equipped with a luer slip connector, and a compartmentalized case.

Connecting Gauge to Cuff

To connect the cuff to the gauge insert the luer slip connector at the end of the bladder tube into the receptacle on the top of the gauge with a firm twist until secure. To remove, reverse procedure.

Using the Trigger Valve

Deflation: Gently squeeze the bottom of the trigger to deflate the cuff. Altering the pressure affects the deflation rate.

Rapid Exhaust: To exhaust remaining air after measurement is completed, squeeze trigger until it clicks and locks in place.

Caution: Never force trigger.


Measurement Procedure

1. **Patient Position:** The patient should sit or lie comfortably. The arm should be fully supported on a flat surface at heart level. (If the arm's position varies, or is not level with the heart, measurement values obtained will not be consistent with the patient's true blood pressure.) Observer should view manometer in a direct line and at eye level to avoid "parallax error".

2. **Apply the Cuff:** ADCUFF™ nylon cuffs with a proprietary Size Guide™ marking system, are specially designed to promote the precise, accurate determination of blood pressure. Index and range markings ensure use of the correct cuff size. The artery mark indicates proper cuff positioning.

Place the cuff over the bare upper arm with artery mark positioned directly over the brachial artery. The bottom edge of the cuff should be positioned approximately one inch above the antecubital fold. Wrap the end of the cuff not containing the bladder around arm snugly and smoothly and engage adhesive strips. To verify a correct fit, check that the INDEX line falls between the two RANGE lines.

3. **Inflate the Cuff:** Palpate the radial artery while inflating the cuff. Be sure to inflate cuff quickly by squeezing bulb rapidly. Inflate cuff 20-30mmHg above the point at which the radial pulse disappears.

4. **Position the Stethoscope:** Position the chestpiece in the antecubital space below the cuff, distal to the brachium. Do not place chestpiece underneath the cuff, as this impedes accurate measurement. Use the bell side of a combination stethoscope for clearest detection of the low pitched Korotkoff (pulse) sounds.

5. **Deflate the Cuff:** Deflate the cuff gradually at the AHA recommended rate of 2-3mmHg per second by gently squeezing the trigger valve mechanism at the rear of the gauge.

6. **Measurement:** Record the onset of Korotkoff sounds (Phase 1) as the systolic pressure, and the disappearance of these sounds (Phase 5) as the diastolic pressure. (Some healthcare professionals prefer recording diastolic 1 and diastolic 2. Diastolic 1 occurs at phase 4).

After measurement is completed, release all remaining air in the cuff by squeezing the trigger until it locks in the rapid "exhaust" position. Remove cuff.


Care and Maintenance

Storage: After measurement, store in carrying case. If unit will not be used for an extended period of time, leave the trigger valve in the full exhaust position.

Manometer: Your ADC® DIAGNOSTIX™ brand palm style aneroid gauge #804 requires minimal care and maintenance.

IMPORTANT: The manometer may be cleaned with a soft cloth but should not be dismantled under any circumstances.

Should the indicator needle of the manometer rest outside of this calibration mark, then the manometer must be recalibrated to within +/-3 mmHg when compared to a certified reference standard. No manometers that have their indicator needle resting outside of this mark are acceptable for use.

NOTE: Store gauge with valve in full exhaust position.

Damaged or broken parts will be replaced as needed at a minimum charge. Refer to the warranty for specific details of warranty coverage.

Inflation System

Consult the enclosed manual entitled "ADCUFF™ SPHYGMOMANOMETER ACCESSORIES" for care and maintenance of inflation system components.

Standards:

Complies with EN 12470-3:2000 and ASTM E 1112-00
The manufacturer recommends a calibration check every 2 years.

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